





# NanoSAM Lab

Scanning Auger  
with Ultimate  
Resolution

- 5 nm SAM Resolution
- 3 nm SEM Resolution
- Depth Profiling and Charge Neutralisation
- Additional Techniques: EBSD, SEMPA, FIB, EBL



SEM  
2  $\mu$ m




# Chemical Analysis with Ultimate Resolution

**The NanoSAM Lab is the ultimate tool for the analysis of small structures. Driven by the unique performance of the UHV Gemini electron column, it guarantees unrivalled resolution below 5 nm in Scanning Auger Microscopy (SAM) and better than 3 nm in SEM.**

In contrast to other Auger tools, the extremely good resolution is not only available at standard 20 keV beam energy, but even at 5 keV the SAM resolution remains below 10 nm. This allows operation in a parameter range where the Auger cross sections are high, and well documented for quantitative analysis.

The NanoSAM Lab is tailored for the fast and efficient acquisition of data on routine samples, while maintaining the flexibility to operate in untypical parameter ranges on challenging materials. Furthermore, the NanoSAM Lab may be extended with additional techniques for a complementary analysis of the key sample characteristics, for example the crystal structure (EBSD) or the magnetic domain structure (SEMPA).



Coloured SEM images of ultra-thin SiN "sheets" on a TEM sample holder, simultaneously acquired with the UHV Gemini Column with in-lens and external secondary electron detector (left and right). Sample by courtesy of: Dr. T. Sekiguchi et al., Nanomaterials Laboratory, NIMS, Japan.





---

## UHV Gemini Column

---

The patented UHV Gemini Column is the electron source with the highest resolution available for UHV operation. Designed for true UHV  $<1 \times 10^{-8}$  Pa, the column is fully bakeable to 180°C and operates without any measurable outgassing: a must for contamination-free research on many materials.

The design of the UHV Gemini Column is the result of a collaboration between Omicron and ZEISS NTS, employing the same electron optics as used in the ZEISS Gemini SEM product range. Both the leading edge SEM software SmartSEM™ and the Gemini control electronics are standard products from Zeiss NTS, incorporated also in the UHV version.

The Auger control software is integrated with the SmartSEM™ software. For efficient data processing an industry standard software is available. Data can also be saved in standard formats such as VAMAS or ASCII.



NanoSAM Lab with sample preparation chamber.

---

## Charge Neutralisation and Depth Profiling

---


Analysis of insulating materials is often impossible with AES, or at least a challenging task. The FIG 05 sputter source enables the operator to compensate negative charges induced by the electron beam with a flood of low-energetic positive ions. This allows analysis of samples which could otherwise not be imaged under stable conditions.

The FIG 05 may also be used for depth profiling with a beam energy ranging from a few 100 eV on delicate samples to 5 keV with high current density for rapid profiling.

## Optional Techniques


- Combined XPS and SAM for in-situ analysis of the elemental composition and chemical functionality.
- Electron BackScatter Diffraction (EBSD) for in-situ analysis of the crystal structure alongside with the elemental composition.
- Focussed Ion Beam (FIB) for tailoring of the sample structure (UHV-Crossbeam with < 10 nm spot size).
- SEM with Polarisation Analysis (SEMPA) for imaging of the magnetic domain structure.
- Electron Beam Lithography (EBL) under cleanest conditions (UHV <  $1 \cdot 10^{-9}$  Pa) for basic research on new resists materials and smallest structures (also for Electron Beam Induced Deposition (EBID)).

For information about the combination of SEM/SAM with SPM and 4-point probes please request our MULTISCAN Lab and UHV NANOPROBE brochures.




**SEMPA:** Image of the magnetic domains of an iron whisker with a typical "fir tree" structure. Each colour indicates a different direction of the magnetisation (4 different domain orientations).


Left: SEMPA detector.



Points, lines and rectangles written with the FIB into a GaAs/AlGaAs heterostructure. Ion-beam-induced secondary electron image.



**EBL:** Test pattern written with the UHV Gemini Column and the Raith ELPHY software.



**FIB:** UHV-compatible 30 kV Ga<sup>+</sup> ion beam source.







Compositional image showing the distribution of silver (yellow), oxygen (red), and silicon (blue).


## NanoSAM Lab at a Glance:

- **UHV Gemini Electron Column:**
  - highest spatial resolution
  - smallest e-beam spot size
  - unique low-energy performance
  - efficient SEM with in-lens SED
- **NanoSAM Electron Analyser:**
  - variable energy resolution
  - excellent sensitivity, multichannel detection
- **Sample Handling:**
  - flexible sample size
  - variable sample temperature (50-500K)
  - sample tilt  $\pm 60^\circ$





Ta	(11.2 nm)
NiFe	(6.8 nm)
AlO <sub>x</sub>	(0.3 nm)
CoFe	(6.8 nm)
Ta	(6.8 nm)
SiO <sub>2</sub>	


Sputter depth profiling experiment on a magnetic tunnel junction. The schematic shows the theoretical structure with 2 magnetic layers (NiFe and CoFe) separated by a thin insulating AlO<sub>x</sub> film. Note that the layer thickness is in the range of the electron mean free path. Nevertheless the different layers are clearly visible, including the extremely thin AlO<sub>x</sub> monolayer.



Silver islands on silicon, Ag MNN, 352 eV.



Gold islands on carbon.



# How to contact us

Headquarters: Omicron NanoTechnology GmbH  
 Limburger Str. 75 • D-65232 Taunusstein • Germany  
 Tel. +49 (0) 61 28 / 987 - 0 • Fax +49 (0) 61 28 / 987 - 185  
 Web: [www.omicron.de](http://www.omicron.de) or [www.omicron-instruments.com](http://www.omicron-instruments.com)  
 e-mail: [info@omicron.de](mailto:info@omicron.de)

## EUROPE

### UK:

Omicron NanoTechnology Ltd.  
 Tel. 01342 33 1000 • Email: [info@omicron.co.uk](mailto:info@omicron.co.uk)

### France:

Omicron NanoTechnology • Tel. 0442 50 68 64  
 Email: [omicron.nanotechnology@wanadoo.fr](mailto:omicron.nanotechnology@wanadoo.fr)

### Italy:

Omicron NanoTechnology GmbH, Roma  
 Tel. 06 35 45 8553 • Email: [omicron.italia@tiscali.it](mailto:omicron.italia@tiscali.it)

### Sweden:

Rowaco AB  
 Tel. 013 138 010 • Email: [kristian@rowaco.se](mailto:kristian@rowaco.se)

### Spain:

TECNOVAC, S.L.  
 Tel. 091 804 1134 • Email: [tecnovac@tecnovac.es](mailto:tecnovac@tecnovac.es)

## AMERICA

### North American Offices:

Omicron NanoTechnology USA  
 Tel. 952 345 5240

Email: [info@omicronUS.com](mailto:info@omicronUS.com)

### North American Offices (North East Regional Office):

Omicron NanoTechnology USA  
 Tel. 724 942 9821

Email: [f.henn@omicronUS.com](mailto:f.henn@omicronUS.com)

### North American Offices (West Regional Office):

Omicron NanoTechnology USA  
 Tel. 303 893 2388

Email: [b.gerace@omicronUS.com](mailto:b.gerace@omicronUS.com)

## ASIA & AUSTRALIA

### Japan:

ULVAC-PHI, Inc.  
 Tel. 0467 85 6522 • Email: [salesjp@phi.com](mailto:salesjp@phi.com)

### South Korea:

Woosin CryoVac Ltd. • Tel. 031 457 5615  
 Email: [woosin@woosinent.co.kr](mailto:woosin@woosinent.co.kr)

### Taiwan:

Omega Scientific Taiwan Ltd. • Tel. 02 8780 5228  
 Email: [omega001@ms3.hinet.net](mailto:omega001@ms3.hinet.net)

### Singapore:

Research Instruments Pte. Ltd.  
 Tel. 6775 7284 • Email: [info@ri.com.sg](mailto:info@ri.com.sg)

### India:


Mack International  
 Tel. 022 2285 5261 • Email: [macktran@vsnl.com](mailto:macktran@vsnl.com)

### China:

Omicron China Office • Tel. 010 8286 7920  
 Email: [germtech@public3.bta.net.cn](mailto:germtech@public3.bta.net.cn)

### Australia:

Thomson Scientific Instruments Pty Ltd.  
 Tel. 03 9663 2738 • Email: [info@tsi.com.au](mailto:info@tsi.com.au)



Sample plate with four contacts.

4-axis sample stage with temperature range <math>< 50\text{ K}</math> up to \pm 60^\circ. Sample dimensions typically  $10 \times 10 \times 5\text{ mm}^3$ , larger on request.